

FOXC2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00470-1

Capture Antibody Information

Catalog Number:
83476-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240403C6
Reactivity:
Human
GenBank:
NM_005251
Immunogen Catalog Number:
Ag33121

Conjugate:
Unconjugated
Full name:
forkhead box C2 (MFH-1, mesenchyme forkhead 1)
Gene ID:
2303

Detection Antibody Information

Catalog Number:
83476-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240403B2
Reactivity:
Human
GenBank:
NM_005251
Immunogen Catalog Number:
Ag33121

Conjugate:
Unconjugated
Full name:
forkhead box C2 (MFH-1, mesenchyme forkhead 1)
Gene ID:
2303

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00470-1 targets FOXC2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FOXC2 Recombinant antibody, PBS Only (Capture) 83476-3-PBS (240403C6). 100 µg. Concentration 1 mg/ml.

Detection antibody: FOXC2 Recombinant antibody, PBS Only (Capture/Detector) 83476-1-PBS (240403B2). 100 µg. Concentration 1 mg/ml.

Alternative FOXC2 matched antibody pairs: MP00470-2, MP00470-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00470-1, FOXC2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83476-3-PBS. Detection antibody: 83476-1-PBS. Standard: Ag33121. Range: 0.313-40 ng/mL