

# HMGB2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00480-1

## Capture Antibody Information

<b>Catalog Number:</b> 83482-2-PBS	<b>Clone ID:</b> 240453G9	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> Human	<b>Full name:</b> high-mobility group box 2
<b>Isotype:</b> IgG	<b>GenBank:</b> BC001063	<b>Gene ID:</b> 3148
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag6135	

## Detection Antibody Information

<b>Catalog Number:</b> 83482-1-PBS	<b>Clone ID:</b> 240453E11	<b>Conjugate:</b> Unconjugated
<b>Host:</b> Rabbit	<b>Reactivity:</b> Human	<b>Full name:</b> high-mobility group box 2
<b>Isotype:</b> IgG	<b>GenBank:</b> BC001063	<b>Gene ID:</b> 3148
<b>Purification Method:</b> Protein A purification	<b>Immunogen Catalog Number:</b> Ag6135	

## Applications

<b>Tested Applications:</b> Cytometric bead array	<b>Range:</b> 3.125-100 ng/mL (Cytometric Bead Array)	<b>Recommended Dilutions:</b> It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

## Product Information

MP00480-1 targets HMGB2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HMGB2 Recombinant antibody, PBS Only (Capture) 83482-2-PBS (240453G9). 100 µg. Concentration 1 mg/ml.

Detection antibody: HMGB2 Recombinant antibody, PBS Only (Detector) 83482-1-PBS (240453E11). 100 µg. Concentration 1 mg/ml.

Alternative HMGB2 matched antibody pairs: MP00480-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP00480-1, HMGB2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83482-2-PBS. Detection antibody: 83482-1-PBS. Standard: Ag6135. Range: 3.125-100 ng/mL.