

HEY2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00482-1

Capture Antibody Information

Catalog Number:
83501-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240474C3
Reactivity:
Human
GenBank:
BC007707
Immunogen Catalog Number:
Ag30988

Conjugate:
Unconjugated
Full name:
hairy/enhancer-of-split related with YRPW motif 2
Gene ID:
23493

Detection Antibody Information

Catalog Number:
83501-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240474F9
Reactivity:
Human
GenBank:
BC007707
Immunogen Catalog Number:
Ag30988

Conjugate:
Unconjugated
Full name:
hairy/enhancer-of-split related with YRPW motif 2
Gene ID:
23493

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00482-1 targets HEY2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HEY2 Recombinant antibody, PBS Only (Capture) 83501-3-PBS (240474C3). 100 µg. Concentration 1 mg/ml.

Detection antibody: HEY2 Recombinant antibody, PBS Only (Detector) 83501-5-PBS (240474F9). 100 µg. Concentration 1 mg/ml.

Alternative HEY2 matched antibody pairs: MP00482-2, MP00482-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00482-1, HEY2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83501-3-PBS. Detection antibody: 83501-5-PBS. Standard: Ag30988. Range: 0.625-40 ng/mL