

For Research Use Only



# HEY2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00482-3

## Capture Antibody Information

**Catalog Number:**  
83501-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240474A6  
**Reactivity:**  
Human  
**GenBank:**  
BC007707  
**Immunogen Catalog Number:**  
Ag30988

**Conjugate:**  
Unconjugated  
**Full name:**  
hairy/enhancer-of-split related with YRPW motif 2  
**Gene ID:**  
23493

## Detection Antibody Information

**Catalog Number:**  
83501-6-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240474F12  
**Reactivity:**  
Human  
**GenBank:**  
BC007707  
**Immunogen Catalog Number:**  
Ag30988

**Conjugate:**  
Unconjugated  
**Full name:**  
hairy/enhancer-of-split related with YRPW motif 2  
**Gene ID:**  
23493

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
1.25-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00482-3 targets HEY2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HEY2 Recombinant antibody, PBS Only (Capture) 83501-1-PBS (240474A6). 100 µg. Concentration 1 mg/ml.

Detection antibody: HEY2 Recombinant antibody, PBS Only (Detector) 83501-6-PBS (240474F12). 100 µg. Concentration 1 mg/ml.

Alternative HEY2 matched antibody pairs: MP00482-1, MP00482-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00482-3, HEV2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83501-1-PBS. Detection antibody: 83501-6-PBS. Standard: Ag30988. Range: 1.25-40 ng/mL