

HTRA2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00491-2

Capture Antibody Information

Catalog Number:
83492-5-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240496G10

Reactivity:
Human

GenBank:
BC000096

Immunogen Catalog Number:
Ag8455

Conjugate:
Unconjugated

Full name:
HtrA serine peptidase 2

Gene ID:
27429

Detection Antibody Information

Catalog Number:
83492-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240496A4

Reactivity:
Human

GenBank:
BC000096

Immunogen Catalog Number:
Ag8455

Conjugate:
Unconjugated

Full name:
HtrA serine peptidase 2

Gene ID:
27429

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-25 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00491-2 targets HTRA2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: HTRA2 Recombinant antibody, PBS Only (Capture) 83492-5-PBS (240496G10). 100 μ g. Concentration 1 mg/mL.

Detection antibody: HTRA2 Recombinant antibody, PBS Only (Detector) 83492-1-PBS (240496A4). 100 μ g. Concentration 1 mg/mL.

Alternative HTRA2 matched antibody pairs: MP00491-1, MP00491-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00491-2, HTRA2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83492-5-PBS. Detection antibody: 83492-1-PBS. Standard: Ag8455. Range: 0.781-25 ng/mL