

For Research Use Only

Mouse IFN-gamma Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00498-1

Capture Antibody Information

Catalog Number:
82768-6-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240262C9

Reactivity:
mouse

GenBank:
BC119063

Conjugate:
Unconjugated

Full name:
interferon gamma

Gene ID:
15978

Detection Antibody Information

Catalog Number:
82768-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240262B6

Reactivity:
mouse

GenBank:
BC119063

Conjugate:
Unconjugated

Full name:
interferon gamma

Gene ID:
15978

Applications

Tested Applications:
Cytometric bead array

Range:
0.78-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00498-1 targets IFN-gamma in immunoassays as a matched antibody pair. Validated in Cytometric bead array. Capture antibody: Mouse IFN gamma Recombinant antibody, PBS Only (Capture) 82768-6-PBS (240262C9). 100 μ g. Concentration 1 mg/ml.

Detection antibody: Mouse IFN gamma Recombinant antibody, PBS Only (Detector) 82768-4-PBS (240262B6). 100 μ g. Concentration 1 mg/ml.

Alternative IFN-gamma matched antibody pairs: MP00498-2, MP01131-1, MP01131-2, MP50092-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00498-1, MOUSE IFN gamma/Ifng Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82768-6-PBS. Detection antibody: 82768-4-PBS. Standard: Eg0334. Range: 0.78-100 ng/mL