

For Research Use Only

# MUC17 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00499-2

## Capture Antibody Information

**Catalog Number:**  
83528-1-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
240351D9

**Reactivity:**  
Human

**GenBank:**  
BC126315

**Conjugate:**  
Unconjugated

**Full name:**  
mucin 17, cell surface associated

**Gene ID:**  
140453

## Detection Antibody Information

**Catalog Number:**  
83528-3-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
240351E2

**Reactivity:**  
Human

**GenBank:**  
BC126315

**Conjugate:**  
Unconjugated

**Full name:**  
mucin 17, cell surface associated

**Gene ID:**  
140453

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.156-20 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00499-2 targets MUC17 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MUC17 Recombinant antibody, PBS Only (Capture) 83528-1-PBS (240351D9). 100 µg. Concentration 1 mg/ml.

Detection antibody: MUC17 Recombinant antibody, PBS Only (Detector) 83528-3-PBS (240351E2). 100 µg. Concentration 1 mg/ml.

Alternative MUC17 matched antibody pairs: MP00499-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00499-2, MUC17 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83528-1-PBS. Detection antibody: 83528-3-PBS. Standard: Eg0413. Range: 0.156-20 ng/mL