

# Integrin beta-5 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00512-3

## Capture Antibody Information

**Catalog Number:**  
83514-4-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240401G1  
**Reactivity:**  
Human  
**GenBank:**  
BC006541  
**Immunogen Catalog Number:**  
Ag29783

**Conjugate:**  
Unconjugated  
**Full name:**  
integrin, beta 5  
**Gene ID:**  
3693

## Detection Antibody Information

**Catalog Number:**  
83514-1-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240401C12  
**Reactivity:**  
Human  
**GenBank:**  
BC006541  
**Immunogen Catalog Number:**  
Ag29783

**Conjugate:**  
Unconjugated  
**Full name:**  
integrin, beta 5  
**Gene ID:**  
3693

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.313-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00512-3 targets Integrin beta-5 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Integrin beta-5 Recombinant antibody, PBS Only (Capture) 83514-4-PBS (240401G1). 100 µg. Concentration 1 mg/ml.

Detection antibody: Integrin beta-5 Recombinant antibody, PBS Only (Detector) 83514-1-PBS (240401C12). 100 µg. Concentration 1 mg/ml.

Alternative Integrin beta-5 matched antibody pairs: MP00512-1, MP00512-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP00512-3, Integrin beta-5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83514-4-PBS. Detection antibody: 83514-1-PBS. Standard: Ag29783. Range: 0.313-40 ng/mL.