

PSMD1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00519-1

Capture Antibody Information

Catalog Number:
83521-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240458A1
Reactivity:
Human
GenBank:
BC094720
Immunogen Catalog Number:
Ag12009

Conjugate:
Unconjugated
Full name:
proteasome (prosome, macropain)
26S subunit, non-ATPase, 1
Gene ID:
5707

Detection Antibody Information

Catalog Number:
83521-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240458B3
Reactivity:
Human
GenBank:
BC094720
Immunogen Catalog Number:
Ag12009

Conjugate:
Unconjugated
Full name:
proteasome (prosome, macropain)
26S subunit, non-ATPase, 1
Gene ID:
5707

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00519-1 targets PSMD1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PSMD1 Recombinant antibody, PBS Only (Capture) 83521-1-PBS (240458A1). 100 µg. Concentration 1 mg/ml.

Detection antibody: PSMD1 Recombinant antibody, PBS Only (Detector) 83521-2-PBS (240458B3). 100 µg. Concentration 1 mg/ml.

Alternative PSMD1 matched antibody pairs: MP00519-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00519-1, PSMD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83521-1-PBS. Detection antibody: 83521-2-PBS. Standard: Ag12009. Range: 0.313-40 ng/mL