

For Research Use Only

MAGI2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00529-1

Capture Antibody Information

Catalog Number:
83561-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240512E9
Reactivity:
Human
GenBank:
BC150277
Immunogen Catalog Number:
Ag18279

Conjugate:
Unconjugated
Full name:
membrane associated guanylate kinase, WW and PDZ domain containing 2
Gene ID:
9863

Detection Antibody Information

Catalog Number:
83561-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240512H3
Reactivity:
Human
GenBank:
BC150277
Immunogen Catalog Number:
Ag18279

Conjugate:
Unconjugated
Full name:
membrane associated guanylate kinase, WW and PDZ domain containing 2
Gene ID:
9863

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00529-1 targets MAGI2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MAGI2 Recombinant antibody, PBS Only (Capture) 83561-3-PBS (240512E9). 100 µg. Concentration 1 mg/mL.

Detection antibody: MAGI2 Recombinant antibody, PBS Only (Detector) 83561-4-PBS (240512H3). 100 µg. Concentration 1 mg/mL.

Alternative MAGI2 matched antibody pairs: MP00529-2, MP00529-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00529-1, MAGI2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83561-3-PBS. Detection antibody: 83561-4-PBS. Standard: Ag18279. Range: 0.781-100 ng/mL.