

DEM1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00536-1

Capture Antibody Information

Catalog Number:
83553-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240504A9
Reactivity:
Human
GenBank:
BC021969
Immunogen Catalog Number:
Ag22666

Conjugate:
Unconjugated
Full name:
defects in morphology 1 homolog (S. cerevisiae)
Gene ID:
64789

Detection Antibody Information

Catalog Number:
83553-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240504A4
Reactivity:
Human
GenBank:
BC021969
Immunogen Catalog Number:
Ag22666

Conjugate:
Unconjugated
Full name:
defects in morphology 1 homolog (S. cerevisiae)
Gene ID:
64789

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00536-1 targets DEM1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DEM1 Recombinant antibody, PBS Only (Capture) 83553-3-PBS (240504A9). 100 µg. Concentration 1 mg/ml.

Detection antibody: DEM1 Recombinant antibody, PBS Only (Detector) 83553-1-PBS (240504A4). 100 µg. Concentration 1 mg/ml.

Alternative DEM1 matched antibody pairs: MP00536-2, MP00536-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00536-1, DEM1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83553-3-PBS. Detection antibody: 83553-1-PBS. Standard: Ag22666. Range: 0.625-40 ng/mL