

For Research Use Only

# SYTL4 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00538-1

## Capture Antibody Information

Catalog Number:  
83564-3-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
240556D9

Reactivity:  
Human

GenBank:  
BC014913

Immunogen Catalog Number:  
Ag2772

Conjugate:  
Unconjugated

Full name:  
synaptotagmin-like 4

Gene ID:  
94121

## Detection Antibody Information

Catalog Number:  
83564-4-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
240556F9

Reactivity:  
Human

GenBank:  
BC014913

Immunogen Catalog Number:  
Ag2772

Conjugate:  
Unconjugated

Full name:  
synaptotagmin-like 4

Gene ID:  
94121

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00538-1 targets SYTL4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SYTL4 Recombinant antibody, PBS Only (Capture) 83564-3-PBS (240556D9). 100 µg. Concentration 1 mg/ml.

Detection antibody: SYTL4 Recombinant antibody, PBS Only (Capture/Detector) 83564-4-PBS (240556F9). 100 µg. Concentration 1 mg/ml.

Alternative SYTL4 matched antibody pairs: MP00538-2, MP00538-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00538-1, SYTL4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83564-3-PBS. Detection antibody: 83564-4-PBS. Standard: Ag2772. Range: 0.625-80 ng/mL