

For Research Use Only

Mouse ICAM-1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00543-3

Capture Antibody Information

Catalog Number:
82827-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240249A10

Reactivity:
Mouse

GenBank:
NM_010493

Conjugate:
Unconjugated

Full name:
intercellular adhesion molecule 1

Gene ID:
15894

Detection Antibody Information

Catalog Number:
82827-4-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240249G1

Reactivity:
Mouse

GenBank:
NM_010493

Conjugate:
Unconjugated

Full name:
intercellular adhesion molecule 1

Gene ID:
15894

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00543-3 targets ICAM-1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse ICAM1 Recombinant antibody, PBS Only (Capture) 82827-2-PBS (240249A10). 100 µg. Concentration 1 mg/ml.

Detection antibody: Mouse ICAM1 Recombinant antibody, PBS Only (Detector) 82827-4-PBS (240249G1). 100 µg. Concentration 1 mg/ml.

Alternative ICAM-1 matched antibody pairs: MP00028-1, MP00028-2, MP00028-3, MP00028-4, MP00028-5, MP00543-1, MP00543-2, MP50085-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00543-3, MOUSE ICAM-1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82827-2-PBS. Detection antibody: 82827-4-PBS. Standard: Eg0482. Range: 0.313-40 ng/mL.