

For Research Use Only

GAPDH Recombinant Matched Antibody Pair, PBS Only



Catalog Number: MP00549-1

Capture Antibody Information

Catalog Number:
81640-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
2L15
Reactivity:
Human
GenBank:
BC004109
Immunogen Catalog Number:
Ag0766

Conjugate:
Unconjugated
Full name:
glyceraldehyde-3-phosphate dehydrogenase
Gene ID:
2597

Detection Antibody Information

Catalog Number:
81640-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
1H18
Reactivity:
Human, mouse, rat
GenBank:
BC004109
Immunogen Catalog Number:
Ag0766

Conjugate:
Unconjugated
Full name:
glyceraldehyde-3-phosphate dehydrogenase
Gene ID:
2597

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00549-1 targets GAPDH in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GAPDH Recombinant antibody, PBS Only (Capture) 81640-3-PBS (2L15). 100 μ g. Concentration 1 mg/mL.

Detection antibody: GAPDH Recombinant antibody, PBS Only (Detector) 81640-5-PBS (1H18). 100 μ g. Concentration 1 mg/mL.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only (BSA and azide free) storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00549-1, GAPDH Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 81640-3-PBS. Detection antibody: 81640-5-PBS. Standard: ag0766. Range: 0.098-100 ng/mL.