

# SEMA3A Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00552-1

## Capture Antibody Information

**Catalog Number:**  
83579-1-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
240509D8

**Reactivity:**  
Human

**GenBank:**  
BC111416

**Immunogen Catalog Number:**  
Ag27173

**Conjugate:**  
Unconjugated

**Full name:**  
sema domain, immunoglobulin domain (Ig), short basic domain, secreted, (semaphorin) 3A

**Gene ID:**  
10371

## Detection Antibody Information

**Catalog Number:**  
83579-2-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
240509G1

**Reactivity:**  
Human

**GenBank:**  
BC111416

**Immunogen Catalog Number:**  
Ag27173

**Conjugate:**  
Unconjugated

**Full name:**  
sema domain, immunoglobulin domain (Ig), short basic domain, secreted, (semaphorin) 3A

**Gene ID:**  
10371

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.313-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00552-1 targets SEMA3A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SEMA3A Recombinant antibody, PBS Only (Capture) 83579-1-PBS (240509D8). 100 µg. Concentration 1 mg/ml.

Detection antibody: SEMA3A Recombinant antibody, PBS Only (Detector) 83579-2-PBS (240509G1). 100 µg. Concentration 1 mg/ml.

Alternative SEMA3A matched antibody pairs: MP00552-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP00552-1, SEMA3A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83579-1-PBS. Detection antibody: 83579-2-PBS. Standard: Ag27173. Range: 0.313-40 ng/mL.