

For Research Use Only

CCDC94 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00553-2

Capture Antibody Information

Catalog Number:
83580-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240517A1

Reactivity:
Human

GenBank:
BC000561

Immunogen Catalog Number:
Ag13764

Conjugate:
Unconjugated

Full name:
coiled-coil domain containing 94

Gene ID:
55702

Detection Antibody Information

Catalog Number:
83580-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240517F11

Reactivity:
Human

GenBank:
BC000561

Immunogen Catalog Number:
Ag13764

Conjugate:
Unconjugated

Full name:
coiled-coil domain containing 94

Gene ID:
55702

Applications

Tested Applications:
Cytometric bead array

Range:
0.78-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00553-2 targets CCDC94 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CCDC94 Recombinant antibody, PBS Only (Capture) 83580-1-PBS (240517A1). 100 μ g. Concentration 1 mg/ml.

Detection antibody: CCDC94 Recombinant antibody, PBS Only (Detector) 83580-2-PBS (240517F11). 100 μ g. Concentration 1 mg/ml.

Alternative CCDC94 matched antibody pairs: MP00553-1, MP00553-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00553-2, CCDC94 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83580-1-PBS. Detection antibody: 83580-2-PBS. Standard: Ag13764. Range: 0.78-100 ng/mL.