

TRIM15 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00555-1

Capture Antibody Information

Catalog Number:
83582-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240565B12
Reactivity:
Human
GenBank:
BC038585
Immunogen Catalog Number:
Ag4533

Conjugate:
Unconjugated
Full name:
tripartite motif-containing 15
Gene ID:
89870

Detection Antibody Information

Catalog Number:
83582-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240565D1
Reactivity:
Human
GenBank:
BC038585
Immunogen Catalog Number:
Ag4533

Conjugate:
Unconjugated
Full name:
tripartite motif-containing 15
Gene ID:
89870

Applications

Tested Applications:
Cytometric bead array

Range:
0.78-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00555-1 targets TRIM15 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TRIM15 Recombinant antibody, PBS Only (Capture) 83582-1-PBS (240565B12). 100 µg. Concentration 1 mg/ml.

Detection antibody: TRIM15 Recombinant antibody, PBS Only (Detector) 83582-3-PBS (240565D1). 100 µg. Concentration 1 mg/ml.

Alternative TRIM15 matched antibody pairs: MP00555-2, MP00555-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00555-1, TRIM15 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83582-1-PBS. Detection antibody: 83582-3-PBS. Standard: Ag4533. Range: 0.78-100 ng/mL