

DMC1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00558-3

Capture Antibody Information

Catalog Number:
83568-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240463E4
Reactivity:
Human
GenBank:
BC035658
Immunogen Catalog Number:
Ag4653

Conjugate:
Unconjugated
Full name:
DMC1 dosage suppressor of mck1 homolog, meiosis-specific homologous recombination (yeast)
Gene ID:
11144

Detection Antibody Information

Catalog Number:
83568-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240463A7
Reactivity:
Human
GenBank:
BC035658
Immunogen Catalog Number:
Ag4653

Conjugate:
Unconjugated
Full name:
DMC1 dosage suppressor of mck1 homolog, meiosis-specific homologous recombination (yeast)
Gene ID:
11144

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00558-3 targets DMC1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: DMC1 Recombinant antibody, PBS Only (Capture) 83568-5-PBS (240463E4). 100 μ g. Concentration 1 mg/ml.

Detection antibody: DMC1 Recombinant antibody, PBS Only (Detector) 83568-1-PBS (240463A7). 100 μ g. Concentration 1 mg/ml.

Alternative DMC1 matched antibody pairs: MP00558-1, MP00558-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00558-3, DMC1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83568-5-PBS. Detection antibody: 83568-1-PBS. Standard: Ag4653. Range: 0.625-80 ng/mL