

TNFB Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00562-3

Capture Antibody Information

Catalog Number:
83573-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240393C2

Reactivity:
Human

GenBank:
BC034729

Conjugate:
Unconjugated

Full name:
lymphotoxin alpha (TNF superfamily, member 1)

Gene ID:
4049

Detection Antibody Information

Catalog Number:
83573-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240393D4

Reactivity:
Human

GenBank:
BC034729

Conjugate:
Unconjugated

Full name:
lymphotoxin alpha (TNF superfamily, member 1)

Gene ID:
4049

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00562-3 targets TNFB in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TNFB Recombinant antibody, PBS Only (Capture/Detector) 83573-1-PBS (240393C2). 100 µg. Concentration 1 mg/mL.

Detection antibody: TNFB Recombinant antibody, PBS Only (Detector) 83573-2-PBS (240393D4). 100 µg. Concentration 1 mg/mL.

Alternative TNFB matched antibody pairs: MP00562-1, MP00562-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00562-3, TNF beta Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83573-1-PBS. Detection antibody: 83573-2-PBS. Standard: Eg0836. Range: 0.313-40 ng/mL.