

GAR1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00568-1

Capture Antibody Information

Catalog Number:
83577-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240537H1
Reactivity:
Human
GenBank:
BC003413
Immunogen Catalog Number:
Ag2282

Conjugate:
Unconjugated
Full name:
GAR1 ribonucleoprotein homolog (yeast)
Gene ID:
54433

Detection Antibody Information

Catalog Number:
83577-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240537C11
Reactivity:
Human
GenBank:
BC003413
Immunogen Catalog Number:
Ag2282

Conjugate:
Unconjugated
Full name:
GAR1 ribonucleoprotein homolog (yeast)
Gene ID:
54433

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00568-1 targets GAR1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GAR1 Recombinant antibody, PBS Only (Capture) 83577-2-PBS (240537H1). 100 µg. Concentration 1 mg/ml.

Detection antibody: GAR1 Recombinant antibody, PBS Only (Detector) 83577-1-PBS (240537C11). 100 µg. Concentration 1 mg/ml.

Alternative GAR1 matched antibody pairs: MP00568-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00568-1, GAR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83577-2-PBS. Detection antibody: 83577-1-PBS. Standard: Ag2282. Range: 0.781-100 ng/mL.