

For Research Use Only

# TCF7L1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00573-1

## Capture Antibody Information

**Catalog Number:**  
83588-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240544F1  
**Reactivity:**  
Human  
**GenBank:**  
BC058894  
**Immunogen Catalog Number:**  
Ag34253

**Conjugate:**  
Unconjugated  
**Full name:**  
transcription factor 7-like 1 (T-cell specific, HMG-box)  
**Gene ID:**  
83439

## Detection Antibody Information

**Catalog Number:**  
83588-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240544G6  
**Reactivity:**  
Human  
**GenBank:**  
BC058894  
**Immunogen Catalog Number:**  
Ag34253

**Conjugate:**  
Unconjugated  
**Full name:**  
transcription factor 7-like 1 (T-cell specific, HMG-box)  
**Gene ID:**  
83439

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
1.25-80 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00573-1 targets TCF7L1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TCF7L1 Recombinant antibody, PBS Only (Capture) 83588-3-PBS (240544F1). 100 µg. Concentration 1 mg/ml.

Detection antibody: TCF7L1 Recombinant antibody, PBS Only (Detector) 83588-2-PBS (240544G6). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00573-1, TCF7L1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83588-3-PBS. Detection antibody: 83588-2-PBS. Standard: Ag34253. Range: 1.25-80 ng/mL