

BARHL2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00582-2

Capture Antibody Information

Catalog Number: 83596-1-PBS	Clone ID: 240558B5	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: BarH-like homeobox 2
Isotype: IgG	GenBank: BC126439	Gene ID: 343472
Purification Method: Protein A purification	Immunogen Catalog Number: Ag21156	

Detection Antibody Information

Catalog Number: 83596-4-PBS	Clone ID: 240558G8	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: BarH-like homeobox 2
Isotype: IgG	GenBank: BC126439	Gene ID: 343472
Purification Method: Protein A purification	Immunogen Catalog Number: Ag21156	

Applications

Tested Applications: Cytometric bead array	Range: 0.313-40 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00582-2 targets BARHL2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BARHL2 Recombinant antibody, PBS Only (Capture) 83596-1-PBS (240558B5). 100 µg. Concentration 1 mg/ml.

Detection antibody: BARHL2 Recombinant antibody, PBS Only (Detector) 83596-4-PBS (240558G8). 100 µg. Concentration 1 mg/ml.

Alternative BARHL2 matched antibody pairs: MP00582-1, MP00582-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00582-2, BARHL2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83596-1-PBS. Detection antibody: 83596-4-PBS. Standard: Ag21156. Range: 0.313-40 ng/mL