

ZBTB7A Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00594-1

Capture Antibody Information

Catalog Number:
83616-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240555B10
Reactivity:
Human
GenBank:
BC113511

Conjugate:
Unconjugated
Full name:
zinc finger and BTB domain containing 7A
Gene ID:
51341

Detection Antibody Information

Catalog Number:
83616-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240555H5
Reactivity:
Human
GenBank:
BC113511

Conjugate:
Unconjugated
Full name:
zinc finger and BTB domain containing 7A
Gene ID:
51341

Applications

Tested Applications:
Cytometric bead array

Range:
39-5000 pg/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00594-1 targets ZBTB7A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ZBTB7A Recombinant antibody, PBS Only (Capture) 83616-1-PBS (240555B10). 100 µg. Concentration 1 mg/mL.

Detection antibody: ZBTB7A Recombinant antibody, PBS Only (Detector) 83616-4-PBS (240555H5). 100 µg. Concentration 1 mg/mL.

Alternative ZBTB7A matched antibody pairs: MP00594-2, MP00594-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00594-1, ZBTB7A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83616-1-PBS. Detection antibody: 83616-4-PBS. Standard: SY00194. Range: 39-5000 pg/mL.