

For Research Use Only

PSEN1-Specific Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00597-1

Capture Antibody Information

Catalog Number:
83618-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240615G10

Reactivity:
Human

GenBank:
NM_000021

Conjugate:
Unconjugated

Full name:
presenilin 1

Gene ID:
5663

Detection Antibody Information

Catalog Number:
83618-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240615B3

Reactivity:
Human

GenBank:
NM_000021

Conjugate:
Unconjugated

Full name:
presenilin 1

Gene ID:
5663

Applications

Tested Applications:
Cytometric bead array

Range:
39-5000 pg/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00597-1 targets PSEN1-Specific in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PSEN1-Specific Recombinant antibody, PBS Only (Capture) 83618-2-PBS (240615G10). 100 µg. Concentration 1 mg/ml.

Detection antibody: PSEN1-Specific Recombinant antibody, PBS Only (Detector) 83618-1-PBS (240615B3). 100 µg. Concentration 1 mg/ml.

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00597-1, PSEN1-Specific Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83618-2-PBS. Detection antibody: 83618-1-PBS. Standard: SY00068. Range: 39-5000 pg/mL