

STEAP2-Specific Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00600-1

Capture Antibody Information

Catalog Number:
83611-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240611F4
Reactivity:
Human
GenBank:
NM_152999

Conjugate:
Unconjugated
Full name:
six transmembrane epithelial antigen of the prostate 2
Gene ID:
261729

Detection Antibody Information

Catalog Number:
83611-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240611B11
Reactivity:
Human
GenBank:
NM_152999

Conjugate:
Unconjugated
Full name:
six transmembrane epithelial antigen of the prostate 2
Gene ID:
261729

Applications

Tested Applications:
Cytometric bead array

Range:
0.039-10 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00600-1 targets STEAP2-Specific in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: STEAP2-Specific Recombinant antibody, PBS Only (Capture) 83611-4-PBS (240611F4). 100 µg. Concentration 1 mg/ml.

Detection antibody: STEAP2-Specific Recombinant antibody, PBS Only (Detector) 83611-2-PBS (240611B11). 100 µg. Concentration 1 mg/ml.

Alternative STEAP2-Specific matched antibody pairs: MP00600-2, MP00600-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00600-1, STEAP2-Specific Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83611-4-PBS. Detection antibody: 83611-2-PBS. Standard: SY00725. Range: 0.039-10 ng/mL.