

Cryptochrome 1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00603-1

Capture Antibody Information

Catalog Number:
83602-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240520H2
Reactivity:
human
GenBank:
BC030519
Immunogen Catalog Number:
Ag4270

Conjugate:
Unconjugated
Full name:
cryptochrome 1 (photolyase-like)
Gene ID:
1407

Detection Antibody Information

Catalog Number:
83602-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240520B4
Reactivity:
human
GenBank:
BC030519
Immunogen Catalog Number:
Ag4270

Conjugate:
Unconjugated
Full name:
cryptochrome 1 (photolyase-like)
Gene ID:
1407

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00603-1 targets Cryptochrome 1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Cryptochrome 1 Recombinant antibody, PBS Only (Capture) 83602-2-PBS (240520H2). 100 µg. Concentration 1 mg/ml.

Detection antibody: Cryptochrome 1 Recombinant antibody, PBS Only (Capture/Detector) 83602-3-PBS (240520B4). 100 µg. Concentration 1 mg/ml.

Alternative Cryptochrome 1 matched antibody pairs: MP00603-2, MP00603-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00603-1, CR1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83602-2-PBS. Detection antibody: 83602-3-PBS. Standard: Ag4270. Range: 0.781-100 ng/mL.