

For Research Use Only

CDK2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00611-1

Capture Antibody Information

Catalog Number:
83635-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240580B9

Reactivity:
Human

GenBank:
BC003065

Immunogen Catalog Number:
Ag17133

Conjugate:
Unconjugated

Full name:
cyclin-dependent kinase 2

Gene ID:
1017

Detection Antibody Information

Catalog Number:
83635-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240580G11

Reactivity:
Human

GenBank:
BC003065

Immunogen Catalog Number:
Ag17133

Conjugate:
Unconjugated

Full name:
cyclin-dependent kinase 2

Gene ID:
1017

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00611-1 targets CDK2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CDK2 Recombinant antibody, PBS Only (Capture) 83635-3-PBS (240580B9). 100 μ g. Concentration 1 mg/mL.

Detection antibody: CDK2 Recombinant antibody, PBS Only (Capture/Detector) 83635-2-PBS (240580G11). 100 μ g. Concentration 1 mg/mL.

Alternative CDK2 matched antibody pairs: MP00611-2, MP00611-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00611-1, CDK2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83635-3-PBS. Detection antibody: 83635-2-PBS. Standard: Ag17133. Range: 0.156-20 ng/mL