

For Research Use Only

CDK2 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00611-4

Capture Antibody Information

Catalog Number:
83635-6-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240580B5

Reactivity:
human

GenBank:
BC003065

Immunogen Catalog Number:
Ag17133

Conjugate:
Unconjugated

Full name:
cyclin-dependent kinase 2

Gene ID:
1017

Detection Antibody Information

Catalog Number:
83635-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240580A9

Reactivity:
human

GenBank:
BC003065

Immunogen Catalog Number:
Ag17133

Conjugate:
Unconjugated

Full name:
cyclin-dependent kinase 2

Gene ID:
1017

Applications

Tested Applications:
Sandwich ELISA

Range:
62.5-4000 pg/mL (Sandwich ELISA)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00611-4 targets CDK2 in immunoassays as a matched antibody pair. Validated in Sandwich ELISA.

Capture antibody: CDK2 Recombinant antibody, PBS Only (Capture) 83635-6-PBS (240580B5). 100 µg. Concentration 1 mg/ml.

Detection antibody: CDK2 Recombinant antibody, PBS Only (Detector) 83635-1-PBS (240580A9). 100 µg. Concentration 1 mg/ml.

Alternative CDK2 matched antibody pairs: MP00611-1, MP00611-2, MP00611-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00611-4, Human CDK2 Recombinant Matched Antibody Pair - PBS only. 83635-6-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag17133. 83635-1-PBS was HRP conjugated as the detection antibody. Range: 62.5-4000 pg/mL.



The mean human CDK2 concentration was determined to be 27.38 ng/mL in Jurkat cell extract based on a 3.0 mg/mL extract load.