

KCND3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00616-2

Capture Antibody Information

Catalog Number:
83640-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240578A10
Reactivity:
human
GenBank:
BC113477
Immunogen Catalog Number:
Ag22183

Conjugate:
Unconjugated
Full name:
potassium voltage-gated channel, Shal-related subfamily, member 3
Gene ID:
3752

Detection Antibody Information

Catalog Number:
83640-5-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240578F12
Reactivity:
human, rat
GenBank:
BC113477
Immunogen Catalog Number:
Ag22183

Conjugate:
Unconjugated
Full name:
potassium voltage-gated channel, Shal-related subfamily, member 3
Gene ID:
3752

Applications

Tested Applications:
Cytometric bead array

Range:
0.78-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00616-2 targets KCND3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: KCND3 Recombinant antibody, PBS Only (Capture) 83640-1-PBS (240578A10). 100 µg. Concentration 1 mg/ml.

Detection antibody: KCND3 Recombinant antibody, PBS Only (Detector) 83640-5-PBS (240578F12). 100 µg. Concentration 1 mg/ml.

Alternative KCND3 matched antibody pairs: MP00616-1, MP00616-3, MP00616-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00616-2, KCND3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83640-1-PBS. Detection antibody: 83640-5-PBS. Standard: Ag22183. Range: 0.78-100 ng/mL