

For Research Use Only

# Involucrin Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00618-2

## Capture Antibody Information

**Catalog Number:**  
83649-4-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
240566G12

**Reactivity:**  
Human

**GenBank:**  
NM\_005547

**Immunogen Catalog Number:**  
Ag29585

**Conjugate:**  
Unconjugated

**Full name:**  
involucrin

**Gene ID:**  
3713

## Detection Antibody Information

**Catalog Number:**  
83649-3-PBS

**Host:**  
Rabbit

**Isotype:**  
IgG

**Purification Method:**  
Protein A purification

**Clone ID:**  
240566F12

**Reactivity:**  
Human

**GenBank:**  
NM\_005547

**Immunogen Catalog Number:**  
Ag29585

**Conjugate:**  
Unconjugated

**Full name:**  
involucrin

**Gene ID:**  
3713

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.156-20 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00618-2 targets Involucrin in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Involucrin Recombinant antibody, PBS Only (Capture) 83649-4-PBS (240566G12). 100 µg. Concentration 1 mg/ml.

Detection antibody: Involucrin Recombinant antibody, PBS Only (Detector) 83649-3-PBS (240566F12). 100 µg. Concentration 1 mg/ml.

Alternative Involucrin matched antibody pairs: MP00618-1, MP00618-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00618-2, Involucrin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83649-4-PBS. Detection antibody: 83649-3-PBS. Standard: Ag29585. Range: 0.156-20 ng/mL.