

For Research Use Only

SOX4 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00619-2

Capture Antibody Information

Catalog Number:
83660-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240625E3

Reactivity:
Human

GenBank:
BC072668

Immunogen Catalog Number:
Ag26583

Conjugate:
Unconjugated

Full name:
SRY (sex determining region Y)-box 4

Gene ID:
6659

Detection Antibody Information

Catalog Number:
83660-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240625H2

Reactivity:
Human

GenBank:
BC072668

Immunogen Catalog Number:
Ag26583

Conjugate:
Unconjugated

Full name:
SRY (sex determining region Y)-box 4

Gene ID:
6659

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00619-2 targets SOX4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SOX4 Recombinant antibody, PBS Only (Capture) 83660-1-PBS (240625E3). 100 µg. Concentration 1 mg/ml.

Detection antibody: SOX4 Recombinant antibody, PBS Only (Detector) 83660-2-PBS (240625H2). 100 µg. Concentration 1 mg/ml.

Alternative SOX4 matched antibody pairs: MP00619-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00619-2, SOX4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83660-1-PBS. Detection antibody: 83660-2-PBS. Standard: Ag26583. Range: 0.625-80 ng/mL