

WFS1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00629-2

Capture Antibody Information

Catalog Number:
83642-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240538B10
Reactivity:
Human
GenBank:
BC030130
Immunogen Catalog Number:
Ag25735

Conjugate:
Unconjugated
Full name:
Wolfram syndrome 1 (wolframin)
Gene ID:
7466

Detection Antibody Information

Catalog Number:
83642-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240538D10
Reactivity:
Human
GenBank:
BC030130
Immunogen Catalog Number:
Ag25735

Conjugate:
Unconjugated
Full name:
Wolfram syndrome 1 (wolframin)
Gene ID:
7466

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00629-2 targets WFS1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: WFS1 Recombinant antibody, PBS Only (Capture) 83642-2-PBS (240538B10). 100 µg. Concentration 1 mg/ml.

Detection antibody: WFS1 Recombinant antibody, PBS Only (Detector) 83642-4-PBS (240538D10). 100 µg. Concentration 1 mg/ml.

Alternative WFS1 matched antibody pairs: MP00629-1, MP00629-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00629-2, WFS1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83642-2-PBS. Detection antibody: 83642-4-PBS. Standard: Ag25735. Range: 0.313-40 ng/mL