

MOSC2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00646-1

Capture Antibody Information

Catalog Number:
83705-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240747F8
Reactivity:
Human
GenBank:
BC011973
Immunogen Catalog Number:
Ag20694

Conjugate:
Unconjugated
Full name:
MOCO sulphurase C-terminal domain containing 2
Gene ID:
54996

Detection Antibody Information

Catalog Number:
83705-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240747A3
Reactivity:
Human
GenBank:
BC011973
Immunogen Catalog Number:
Ag20694

Conjugate:
Unconjugated
Full name:
MOCO sulphurase C-terminal domain containing 2
Gene ID:
54996

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00646-1 targets MOSC2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MOSC2 Recombinant antibody, PBS Only (Capture) 83705-3-PBS (240747F8). 100 µg. Concentration 1 mg/ml.

Detection antibody: MOSC2 Recombinant antibody, PBS Only (Detector) 83705-1-PBS (240747A3). 100 µg. Concentration 1 mg/ml.

Alternative MOSC2 matched antibody pairs: MP00646-2, MP00646-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00646-1, MOSC2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83705-3-PBS. Detection antibody: 83705-1-PBS. Standard: Ag20694. Range: 0.313-20 ng/mL