

RAB37 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00690-2

Capture Antibody Information

Catalog Number:
83740-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240838A6

Reactivity:
Human

GenBank:
BC016615

Immunogen Catalog Number:
Ag3717

Conjugate:
Unconjugated

Full name:
RAB37, member RAS oncogene family

Gene ID:
326624

Detection Antibody Information

Catalog Number:
83740-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240838D9

Reactivity:
Human

GenBank:
BC016615

Immunogen Catalog Number:
Ag3717

Conjugate:
Unconjugated

Full name:
RAB37, member RAS oncogene family

Gene ID:
326624

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00690-2 targets RAB37 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RAB37 Recombinant antibody, PBS Only (Capture/Detector) 83740-1-PBS (240838A6). 100 µg. Concentration 1 mg/ml.

Detection antibody: RAB37 Recombinant antibody, PBS Only (Detector) 83740-2-PBS (240838D9). 100 µg. Concentration 1 mg/ml.

Alternative RAB37 matched antibody pairs: MP00690-1, MP00690-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00690-2, RAB37 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83740-1-PBS. Detection antibody: 83740-2-PBS. Standard: Ag3717. Range: 0.625-80 ng/mL