

AOX Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00709-1

Capture Antibody Information

Catalog Number:
83731-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240754G9
Reactivity:
Human
GenBank:
BC010425
Immunogen Catalog Number:
Ag1405

Conjugate:
Unconjugated
Full name:
acyl-Coenzyme A oxidase 1, palmitoyl
Gene ID:
51

Detection Antibody Information

Catalog Number:
83731-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240754D4
Reactivity:
Human
GenBank:
BC010425
Immunogen Catalog Number:
Ag1405

Conjugate:
Unconjugated
Full name:
acyl-Coenzyme A oxidase 1, palmitoyl
Gene ID:
51

Applications

Tested Applications:
Cytometric bead array

Range:
0.78-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00709-1 targets AOX in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: AOX Recombinant antibody, PBS Only (Capture) 83731-4-PBS (240754G9). 100 µg. Concentration 1 mg/mL.

Detection antibody: AOX Recombinant antibody, PBS Only (Capture/Detector) 83731-3-PBS (240754D4). 100 µg. Concentration 1 mg/mL.

Alternative AOX matched antibody pairs: MP00709-2, MP00709-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00709-1, AOX Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83731-4-PBS. Detection antibody: 83731-3-PBS. Standard: Ag1405. Range: 0.78-100 ng/mL