

MIB1 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00718-2

Capture Antibody Information

Catalog Number:
83718-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240676E6
Reactivity:
Human
GenBank:
BC022403
Immunogen Catalog Number:
Ag2507

Conjugate:
Unconjugated
Full name:
mindbomb homolog 1 (Drosophila)
Gene ID:
57534

Detection Antibody Information

Catalog Number:
83718-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240676E10
Reactivity:
Human
GenBank:
BC022403
Immunogen Catalog Number:
Ag2507

Conjugate:
Unconjugated
Full name:
mindbomb homolog 1 (Drosophila)
Gene ID:
57534

Applications

Tested Applications:
Cytometric bead array

Range:
1.25-80 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00718-2 targets MIB1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MIB1 Recombinant antibody, PBS Only (Capture) 83718-2-PBS (240676E6). 100 µg. Concentration 1 mg/ml.

Detection antibody: MIB1 Recombinant antibody, PBS Only (Capture/Detector) 83718-4-PBS (240676E10). 100 µg. Concentration 1 mg/ml.

Alternative MIB1 matched antibody pairs: MP00718-1, MP00718-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00718-2, MIB1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83718-2-PBS. Detection antibody: 83718-4-PBS. Standard: Ag2507. Range: 1.25-80 ng/mL