

# FOXK2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00722-1

## Capture Antibody Information

**Catalog Number:**  
83735-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240692E8  
**Reactivity:**  
Human  
**GenBank:**  
NM\_004514.3  
**Immunogen Catalog Number:**  
Ag33625

**Conjugate:**  
Unconjugated  
**Full name:**  
forkhead box K2  
**Gene ID:**  
3607

## Detection Antibody Information

**Catalog Number:**  
83735-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240692B7  
**Reactivity:**  
Human  
**GenBank:**  
NM\_004514.3  
**Immunogen Catalog Number:**  
Ag33625

**Conjugate:**  
Unconjugated  
**Full name:**  
forkhead box K2  
**Gene ID:**  
3607

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.313-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00722-1 targets FOXK2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FOXK2 Recombinant antibody, PBS Only (Capture) 83735-2-PBS (240692E8). 100 µg. Concentration 1 mg/ml.

Detection antibody: FOXK2 Recombinant antibody, PBS Only (Capture/Detector) 83735-3-PBS (240692B7). 100 µg. Concentration 1 mg/ml.

Alternative FOXK2 matched antibody pairs: MP00722-2, MP00722-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**Storage buffer:**  
PBS only

## Selected Validation Data



Cytometric bead array standard curve of MP00722-1, FOXP2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83735-2-PBS. Detection antibody: 83735-3-PBS. Standard: Ag33625. Range: 0.313-40 ng/mL