

FOXK2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number: MP00722-3

Capture Antibody Information

Catalog Number:
83735-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240692E8
Reactivity:
Human
GenBank:
NM_004514.3
Immunogen Catalog Number:
Ag33625

Conjugate:
Unconjugated
Full name:
forkhead box K2
Gene ID:
3607

Detection Antibody Information

Catalog Number:
83735-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240692A7
Reactivity:
Human
GenBank:
NM_004514.3
Immunogen Catalog Number:
Ag33625

Conjugate:
Unconjugated
Full name:
forkhead box K2
Gene ID:
3607

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00722-3 targets FOXK2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FOXK2 Recombinant antibody, PBS Only (Capture) 83735-2-PBS (240692E8). 100 µg. Concentration 1 mg/ml.

Detection antibody: FOXK2 Recombinant antibody, PBS Only (Detector) 83735-1-PBS (240692A7). 100 µg. Concentration 1 mg/ml.

Alternative FOXK2 matched antibody pairs: MP00722-1, MP00722-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00722-3, FOXP2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83735-2-PBS. Detection antibody: 83735-1-PBS. Standard: Ag33625. Range: 0.313-40 ng/mL