

YAF2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00726-1

Capture Antibody Information

Catalog Number:
83798-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240834B9
Reactivity:
Human
GenBank:
BC037777
Immunogen Catalog Number:
Ag5023

Conjugate:
Unconjugated
Full name:
YY1 associated factor 2
Gene ID:
10138

Detection Antibody Information

Catalog Number:
83798-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240834F11
Reactivity:
Human
GenBank:
BC037777
Immunogen Catalog Number:
Ag5023

Conjugate:
Unconjugated
Full name:
YY1 associated factor 2
Gene ID:
10138

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00726-1 targets YAF2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: YAF2 Recombinant antibody, PBS Only (Capture) 83798-3-PBS (240834B9). 100 μ g. Concentration 1 mg/ml.

Detection antibody: YAF2 Recombinant antibody, PBS Only (Detector) 83798-2-PBS (240834F11). 100 μ g. Concentration 1 mg/ml.

Alternative YAF2 matched antibody pairs: MP00726-2, MP00726-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00726-1, YAF2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83798-3-PBS. Detection antibody: 83798-2-PBS. Standard: Ag5023. Range: 0.625-40 ng/mL