

For Research Use Only

Rat RAGE Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00734-3

Capture Antibody Information

Catalog Number:
83759-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240706D1
Reactivity:
Rat
GenBank:

Conjugate:
Unconjugated
Full name:
advanced glycosylation end product-specific receptor
Gene ID:
81722

Detection Antibody Information

Catalog Number:
83759-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240706D12
Reactivity:
Rat
GenBank:

Conjugate:
Unconjugated
Full name:
advanced glycosylation end product-specific receptor
Gene ID:
81722

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00734-3 targets RAGE in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Rat RAGE Recombinant antibody, PBS Only (Capture/Detector) 83759-3-PBS (240706D1). 100 μ g. Concentration 1 mg/ml.

Detection antibody: Rat RAGE Recombinant antibody, PBS Only (Capture/Detector) 83759-2-PBS (240706D12). 100 μ g. Concentration 1 mg/ml.

Alternative RAGE matched antibody pairs: MP00693-1, MP00693-2, MP00693-3, MP00734-1, MP00734-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00734-3, RAT RAGE Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83759-3-PBS. Detection antibody: 83759-2-PBS. Standard: Eg1135. Range: 0.156-20 ng/mL.