

WDR21C Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00738-2

Capture Antibody Information

Catalog Number:
83802-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240402C11
Reactivity:
Human
GenBank:
BC111778
Immunogen Catalog Number:
Ag30748

Conjugate:
Unconjugated
Full name:
WD repeat domain 21C
Gene ID:
138009

Detection Antibody Information

Catalog Number:
83802-4-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240402H3
Reactivity:
Human
GenBank:
BC111778
Immunogen Catalog Number:
Ag30748

Conjugate:
Unconjugated
Full name:
WD repeat domain 21C
Gene ID:
138009

Applications

Tested Applications:
Cytometric bead array

Range:
1.56-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00738-2 targets WDR21C in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: WDR21C Recombinant antibody, PBS Only (Capture) 83802-1-PBS (240402C11). 100 µg. Concentration 1 mg/ml.

Detection antibody: WDR21C Recombinant antibody, PBS Only (Detector) 83802-4-PBS (240402H3). 100 µg. Concentration 1 mg/ml.

Alternative WDR21C matched antibody pairs: MP00738-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00738-2, WDR21C Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83802-1-PBS. Detection antibody: 83802-4-PBS. Standard: Ag30748. Range: 1.56-100 ng/mL.