

BAG3 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00740-2

Capture Antibody Information

Catalog Number:
83779-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240860C9
Reactivity:
Human
GenBank:
BC006418
Immunogen Catalog Number:
Ag0956

Conjugate:
Unconjugated
Full name:
BCL2-associated athanogene 3
Gene ID:
9531

Detection Antibody Information

Catalog Number:
83779-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240860D4
Reactivity:
Human
GenBank:
BC006418
Immunogen Catalog Number:
Ag0956

Conjugate:
Unconjugated
Full name:
BCL2-associated athanogene 3
Gene ID:
9531

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-40 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00740-2 targets BAG3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BAG3 Recombinant antibody, PBS Only (Capture) 83779-1-PBS (240860C9). 100 µg. Concentration 1 mg/ml.

Detection antibody: BAG3 Recombinant antibody, PBS Only (Detector) 83779-3-PBS (240860D4). 100 µg. Concentration 1 mg/ml.

Alternative BAG3 matched antibody pairs: MP00740-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00740-2, BAG3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83779-1-PBS. Detection antibody: 83779-3-PBS. Standard: Ag0956. Range: 0.625-40 ng/mL