

MITF Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00746-1

Capture Antibody Information

Catalog Number:
83803-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240815B4

Reactivity:
human

GenBank:
BC012503

Immunogen Catalog Number:
Ag3679

Conjugate:
Unconjugated

Full name:
microphthalmia-associated
transcription factor

Gene ID:
4286

Detection Antibody Information

Catalog Number:
83803-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240815E7

Reactivity:
human

GenBank:
BC012503

Immunogen Catalog Number:
Ag3679

Conjugate:
Unconjugated

Full name:
microphthalmia-associated
transcription factor

Gene ID:
4286

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00746-1 targets MITF in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: MITF Recombinant antibody, PBS Only (Capture/Detector) 83803-1-PBS (240815B4). 100 µg. Concentration 1 mg/ml.

Detection antibody: MITF Recombinant antibody, PBS Only (Detector) 83803-2-PBS (240815E7). 100 µg. Concentration 1 mg/ml.

Alternative MITF matched antibody pairs: MP00746-2, MP00746-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:


Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00746-1, MITF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83803-1-PBS. Detection antibody: 83803-2-PBS. Standard: Ag3679. Range: 0.781-100 ng/mL.