

For Research Use Only

# ISLR Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00747-2

## Capture Antibody Information

Catalog Number:  
83782-3-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
240727A8

Reactivity:  
Human

GenBank:  
BC022478

Immunogen Catalog Number:  
Ag19446

Conjugate:  
Unconjugated

Full name:  
immunoglobulin superfamily  
containing leucine-rich repeat

Gene ID:  
3671

## Detection Antibody Information

Catalog Number:  
83782-4-PBS

Host:  
Rabbit

Isotype:  
IgG

Purification Method:  
Protein A purification

Clone ID:  
240727C1

Reactivity:  
Human

GenBank:  
BC022478

Immunogen Catalog Number:  
Ag19446

Conjugate:  
Unconjugated

Full name:  
immunoglobulin superfamily  
containing leucine-rich repeat

Gene ID:  
3671

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.313-20 ng/mL (Cytometric Bead  
Array)

Recommended Dilutions:  
It is recommended that this reagent  
should be titrated in each testing  
system to obtain optimal results.

## Product Information

MP00747-2 targets ISLR in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ISLR Recombinant antibody, PBS Only (Capture) 83782-3-PBS (240727A8). 100  $\mu$ g. Concentration 1 mg/mL.

Detection antibody: ISLR Recombinant antibody, PBS Only (Detector) 83782-4-PBS (240727C1). 100  $\mu$ g. Concentration 1 mg/mL.

Alternative ISLR matched antibody pairs: MP00747-1, MP00747-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.  
Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00747-2, ISLR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83782-3-PBS. Detection antibody: 83782-4-PBS. Standard: Ag19446. Range: 0.313-20 ng/mL