

RXRA Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00759-1

Capture Antibody Information

Catalog Number: 83796-3-PBS	Clone ID: 240878B1	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: retinoid X receptor, alpha
Isotype: IgG	GenBank: BC110998	Gene ID: 6256
Purification Method: Protein A purification	Immunogen Catalog Number: Ag15651	

Detection Antibody Information

Catalog Number: 83796-1-PBS	Clone ID: 240878A4	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: retinoid X receptor, alpha
Isotype: IgG	GenBank: BC110998	Gene ID: 6256
Purification Method: Protein A purification	Immunogen Catalog Number: Ag15651	

Applications

Tested Applications: Cytometric bead array	Range: 0.625-80 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP00759-1 targets RXRA in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RXRA Recombinant antibody, PBS Only (Capture) 83796-3-PBS (240878B1). 100 µg. Concentration 1 mg/ml.

Detection antibody: RXRA Recombinant antibody, PBS Only (Detector) 83796-1-PBS (240878A4). 100 µg. Concentration 1 mg/ml.

Alternative RXRA matched antibody pairs: MP00759-2, MP00759-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00759-1, RXRA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83796-3-PBS. Detection antibody: 83796-1-PBS. Standard: Ag15651. Range: 0.625-80 ng/mL