

For Research Use Only

IL33 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00763-1

Capture Antibody Information

Catalog Number:
83765-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240719G4

Reactivity:
Human

GenBank:
BC047085

Conjugate:
Unconjugated

Full name:
interleukin 33

Gene ID:
90865

Detection Antibody Information

Catalog Number:
83765-3-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
240719F5

Reactivity:
Human

GenBank:
BC047085

Conjugate:
Unconjugated

Full name:
interleukin 33

Gene ID:
90865

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00763-1 targets IL33 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IL33 Recombinant antibody, PBS Only (Capture) 83765-2-PBS (240719G4). 100 μ g. Concentration 1 mg/ml.

Detection antibody: IL33 Recombinant antibody, PBS Only (Detector) 83765-3-PBS (240719F5). 100 μ g. Concentration 1 mg/ml.

Alternative IL33 matched antibody pairs: MP00763-2, MP00763-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00763-1, IL-33 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83765-2-PBS. Detection antibody: 83765-3-PBS. Standard: Eg0313. Range: 1.563-100 ng/mL.