

For Research Use Only

Mouse Cdkn2a Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00764-2

Capture Antibody Information

Catalog Number:
83797-2-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240935F9
Reactivity:
mouse
GenBank:
NM_001040654
Immunogen Catalog Number:
Ag28927

Conjugate:
Unconjugated
Full name:
cyclin-dependent kinase inhibitor 2A
Gene ID:
12578

Detection Antibody Information

Catalog Number:
83797-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240935E2
Reactivity:
mouse
GenBank:
NM_001040654
Immunogen Catalog Number:
Ag28927

Conjugate:
Unconjugated
Full name:
cyclin-dependent kinase inhibitor 2A
Gene ID:
12578

Applications

Tested Applications:
Cytometric bead array

Range:
3.13-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00764-2 targets Cdkn2a in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Mouse Cdkn2a/p16INK4a Recombinant antibody, PBS Only (Capture) 83797-2-PBS (240935F9). 100 µg. Concentration 1 mg/ml.

Detection antibody: Mouse Cdkn2a/p16INK4a Recombinant antibody, PBS Only (Detector) 83797-1-PBS (240935E2). 100 µg. Concentration 1 mg/ml.

Alternative Cdkn2a matched antibody pairs: MP00764-1, MP00764-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00764-2, MOUSE Cdkn2a Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83797-2-PBS. Detection antibody: 83797-1-PBS. Standard: Ag28927. Range: 3.13-100 ng/mL.