

For Research Use Only

# SETDB1 Recombinant Matched Antibody Pair, PBS Only



Catalog Number:MP00766-2

## Capture Antibody Information

**Catalog Number:**  
83772-2-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240805F5  
**Reactivity:**  
human  
**GenBank:**  
BC009362  
**Immunogen Catalog Number:**  
Ag1725

**Conjugate:**  
Unconjugated  
**Full name:**  
SET domain, bifurcated 1  
**Gene ID:**  
9869

## Detection Antibody Information

**Catalog Number:**  
83772-3-PBS  
**Host:**  
Rabbit  
**Isotype:**  
IgG  
**Purification Method:**  
Protein A purification

**Clone ID:**  
240805E8  
**Reactivity:**  
human  
**GenBank:**  
BC009362  
**Immunogen Catalog Number:**  
Ag1725

**Conjugate:**  
Unconjugated  
**Full name:**  
SET domain, bifurcated 1  
**Gene ID:**  
9869

## Applications

**Tested Applications:**  
Cytometric bead array

**Range:**  
0.625-40 ng/mL (Cytometric Bead Array)

**Recommended Dilutions:**  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP00766-2 targets SETDB1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SETDB1 Recombinant antibody, PBS Only (Capture) 83772-2-PBS (240805F5). 100 µg. Concentration 1 mg/ml.

Detection antibody: SETDB1 Recombinant antibody, PBS Only (Detector) 83772-3-PBS (240805E8). 100 µg. Concentration 1 mg/ml.

Alternative SETDB1 matched antibody pairs: MP00766-1, MP00766-3, MP00766-4

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
**Storage buffer:**  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00766-2, SETDB1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83772-2-PBS. Detection antibody: 83772-3-PBS. Standard: Ag1725. Range: 0.625-40 ng/mL