

SON Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00772-2

Capture Antibody Information

Catalog Number:
83787-3-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240845D7
Reactivity:
Human
GenBank:
NM_001291411
Immunogen Catalog Number:
Ag27620

Conjugate:
Unconjugated
Full name:
SON DNA binding protein
Gene ID:
6651

Detection Antibody Information

Catalog Number:
83787-1-PBS
Host:
Rabbit
Isotype:
IgG
Purification Method:
Protein A purification

Clone ID:
240845C12
Reactivity:
Human
GenBank:
NM_001291411
Immunogen Catalog Number:
Ag27620

Conjugate:
Unconjugated
Full name:
SON DNA binding protein
Gene ID:
6651

Applications

Tested Applications:
Cytometric bead array

Range:
0.313-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP00772-2 targets SON in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SON Recombinant antibody, PBS Only (Capture) 83787-3-PBS (240845D7). 100 μ g. Concentration 1 mg/mL.

Detection antibody: SON Recombinant antibody, PBS Only (Detector) 83787-1-PBS (240845C12). 100 μ g. Concentration 1 mg/mL.

Alternative SON matched antibody pairs: MP00772-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00772-2, SON Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83787-3-PBS. Detection antibody: 83787-1-PBS. Standard: Ag27620. Range: 0.313-20 ng/mL