

TNFR2 Recombinant Matched Antibody Pair, PBS Only

Catalog Number:MP00777-1

Capture Antibody Information

Catalog Number:
83101-2-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230328E4

Reactivity:
human

GenBank:
BC052977

Conjugate:
Unconjugated

Full name:
tumor necrosis factor receptor
superfamily, member 1B

Gene ID:
7133

Detection Antibody Information

Catalog Number:
98010-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
230328B3

Reactivity:
human

GenBank:
BC052977

Conjugate:
Unconjugated

Full name:
tumor necrosis factor receptor
superfamily, member 1B

Gene ID:
7133

Applications

Tested Applications:
Cytometric bead array

Range:
0.625-20 ng/mL (Cytometric Bead
Array)

Recommended Dilutions:
It is recommended that this reagent
should be titrated in each testing
system to obtain optimal results.

Product Information

MP00777-1 targets TNFR2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TNFR2 Recombinant antibody, PBS Only (Capture) 83101-2-PBS (230328E4). 100 µg. Concentration 1 mg/ml.

Detection antibody: Anti-Human TNFR2/CD120b Rabbit Recombinant Antibody, PBS Only 98010-1-PBS (230328B3). 100 µg. Concentration 1 mg/ml.

Alternative TNFR2 matched antibody pairs: MP00777-2, MP00777-3, MP00777-4, MP50139-1

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.


Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP00777-1, TNFRSF1B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 83101-2-PBS. Detection antibody: 98010-1-PBS. Standard: Eg0086. Range: 0.625-20 ng/mL.